

# GOJO® SMARTLINK™ Activity Monitoring System

Community Level Hand Hygiene Compliance System

## Did You Know?

Hand hygiene compliance rates average

**less than 50%** nationally<sup>1</sup>

Initial research indicates our Solution enables

a **92% increase** in hand hygiene compliance rates<sup>2</sup>

A 1.0% increase in hand hygiene compliance resulted in annual savings of

**\$39,650**<sup>3</sup>



Improved patient outcomes for a higher standard of care.

GOJO SMARTLINK Activity Monitoring System can monitor and measure compliance at the community level. It can seamlessly convert real-time hand hygiene opportunity and activity data into actionable, usable information by floor, unit or room.

- Easy-to-install GOJO-licensed hardware and customizable software system
- Efficient 24/7 gathering of objective data
- Delivers real-time actionable data
- Works with GOJO products, dispensing technology and support to build a successful program

For more information, call 800-321-9647, visit [healthcare.GOJO.com/SMARTLINK](http://healthcare.GOJO.com/SMARTLINK) or e-mail [healthcare@gojo.com](mailto:healthcare@gojo.com).



<sup>1</sup>Herbert C, Weber SG. Common approaches to the control of multidrug-resistant organisms other than methicillin-resistant Staphylococcus aureus (MRSA). *Infect Dis Clin North Am.* 2011 Mar;25(1): 181-200. Epub 2010 Dec 17.

<sup>2</sup>GOJO Pilot. Fort Worth, Texas Hospital. Two units studied June – October 2012.

<sup>3</sup>To a 200-bed hospital. Cummings KL, Anderson DJ, Kaye KS. Hand hygiene noncompliance and the cost of hospital-acquired methicillin-resistant Staphylococcus aureus infection. *Infection Control and Hospital Epidemiology*, April 2010.

# GOJO® SMARTLINK™ Activity Monitoring System

## community level system

This system accurately captures and estimates hand hygiene compliance at the community level within your facility. It can be configured to monitor and measure hand hygiene compliance by floor, unit or room.

$$\frac{\text{ACTIVITIES}}{\text{OPPORTUNITIES}} = \text{COMPLIANCE RATE\%}$$



step 1

"Cleaning in" before entering a patient room is recorded as a hand hygiene "activity" and transmitted.



step 2

Entry into the patient room is captured as a hand hygiene "opportunity" and transmitted.



step 3

Patient care is provided.



step 4

"Cleaning out" is recorded as a hand hygiene "activity" and transmitted.



step 5

Exiting the patient room is captured as a hand hygiene "opportunity" and transmitted.



step 6

SMARTLINK software uses real-time "activity" and "opportunity" data to calculate hand hygiene compliance rates and generate reports.



### Here's What's Included

- GOJO SMARTLINK Technology modules for dispensers
- One activity counter installed per room
- Network communication components installed as needed

### Our Software System Offers

- Easy-to-use graphical user interface
- Data reported in real time (within 1 minute) and continually updated
- Data that can be viewed by time and place—even to the room level
- Data outputs that can be customized by opportunities
- Custom slide shows to motivate staff
- Viewing on multiple devices
- Alerts for sensor failures, low batteries or other device issues

### GOJO SMARTLINK Clinician-Based Support (Optional)

- Clinician-directed program that helps support and sustain your hand hygiene compliance improvement program
- Brings to life the practices you need to continue improvement
- Based on Six Sigma improvement processes
- Can be customized to meet the needs of your facility



If you're ready to discuss your hand hygiene compliance goals, contact your GOJO sales representative today.  
**Visit our website at [healthcare.GOJO.com/SMARTLINK](http://healthcare.GOJO.com/SMARTLINK).**

#### Information for Your IT Department:

- Transmits small packets of data (2–4 bits)
- Data is only outbound
- Operates on 900MHz for data transmission
- Does not require extensive network integration
- Battery-powered wireless devices
- No HIPAA data is transmitted

GOJO Industries, Inc.  
One GOJO Plaza, Suite 500  
P.O. Box 991  
Akron, OH 44309-0991  
Tel: 1-330-255-6000  
Toll-free:  
1-800-321-9647  
Fax: 1-800-FAX-GOJO

GOJO Industries - Europe, Ltd.\*  
Units 5 & 6  
Stratus Park  
Brinklow  
Milton Keynes  
MK10 0DE UK  
Tel: +44 (0) 1908 588444  
Fax: +44 (0) 1908 588445

GOJO France SARL\*  
90 rue de Paris  
59800 LILLE  
France  
Tél: +33 (0)3 20 30 33 38  
Fax: +33 (0)3 20 48 28 71

GOJO América Latina, Ltda.\*  
Av. Nossa Sra. do Bom Sucesso,  
3344  
Condomínio Industrial Turn Key  
Módulo 1  
Pindamonhangaba - SP - Brasil  
CEP: 12.420-010  
Tel/Fax: +55 12 3644 2600

GOJO Japan, Inc.\*  
MIF Building  
3rd Floor  
1-6-6, Uchikanda  
Chiyoda-ku, Tokyo 101-0047  
Japan  
Tel: +81-3-5280-4807  
Fax: +81-3-5280-4843

GOJO Australasia Pty Ltd\*  
Suite 106,  
460 Pacific Highway  
St Leonards, NSW 2065  
Australia  
Tel: +61 2 9016 3882  
Fax: +61 2 9437 5571